

Ver 1.0

g

Handbook for

■ Soil DNA mini

exgene™

DNA PURIFICATION HANDBOOK


GeneAll

Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.) : tech@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

AmpONE™, Exfection™, Exgene™, Expin™, Exprep™, EzClear™, EzSep™, EzPure™, GenEx™, Hybrid-Q™, DirEx™, Allspin™, RiboEx™, Riboclear™, Ribospin™ are trademarks of GeneAll Biotechnology co., ltd.

© 2012 GeneAll Biotechnology, all right reserved.

This protocol handbook is included in :

GeneAll® Exgene™ Soil DNA mini (I14-I50)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

Sample Pulverization step

Add up to 500 mg of soil sample to a Powerbead™ tube.
Add 550 ul of Buffer SL.
Pulverize the sample.
Centrifuge at $\geq 10,000 \times g$ for 10 minutes.

Inhibitor removal step

Transfer the supernatant to a 1.5 ml tube.
Add 50 ul of buffer RH.
Add 300 ul of buffer PD and mix well.
Centrifuge at $\geq 10,000 \times g$ for 5 minutes.

DNA binding step

Transfer the supernatant to a 2 ml tube.
Add 900 ul of buffer TB.
Apply the mixture into a mini spin column and centrifuge at $\geq 10,000 \times g$ for 30 seconds.


Washing step

Add 500 ul of buffer NW and
Centrifuge at $\geq 10,000 \times g$ for 30 seconds.
Centrifuge at $\geq 10,000 \times g$ for 1 minute.

DNA elution

Add ~50 ul of Buffer EB to the center of the membrane.
Centrifuge at $\geq 10,000 \times g$ for 1 minute.

Brief protocol



INDEX

| | |
|------------------------|-----------|
| Index | 05 |
| Kit Contents | 06 |
| Materials not provided | |
| Quality Control | 07 |
| Storage Conditions | |
| Precautions | |
| Product disclaimer | |
| Product specifications | 08 |
| Product description | 09 |
| Protocol | 10 |
| Troubleshooting Guide | 12 |

GeneAll®
Exgene™ Soil DNA mini

KIT CONTENTS

| Components | Quantity | Storage |
|---|-----------------|---------------------|
| Buffer SL | 30 ml | Room temperature |
| Buffer RH | 3 ml | |
| Buffer PD | 17 ml | |
| Buffer TB | 50 ml | |
| Buffer NW | 30 ml | |
| Buffer EB | 15 ml | |
| Powerbead™ tube | 50 | |
| GeneAll® Column type G (with collection tube) | 50 | |
| 1.5 ml tube | 100 | |
| 2.0 ml tube | 50 | |

MATERIALS NOT PROVIDED

Disposable material

- Pipet tips
- Disposable gloves

Equipment

- Precellys®24 (Bertin, France) equipment or any equivalent
- Microcentrifuge
- Suitable protector (ex; lab coat, disposable gloves, goggles, etc)

QUALITY CONTROL

GeneAll® Exgene™ Soil DNA mini is manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. For consistency of product, the quality certification process is carried out from lot to lot thoroughly and only the qualified is approved to be delivered.

STORAGE CONDITIONS

GeneAll® Exgene™ Soil DNA mini should be stored at room temperature (15 ~ 25°C). But prolonged storage at high temperature over 30°C can reduce the performance of the kit.

In cold ambient condition, buffer RH and TB may exhibit salt precipitation and this will cause reduction of DNA recover-yields. If so, heat the bottle with occasional swirling in 37°C water bath until completely dissolved.

All components are stable for 1 year.

Keep out of direct sunlight.

PRECAUTIONS

The buffers included in GeneAll® Exgene™ Soil DNA mini contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions. In case of contact, wash immediately with plenty of water and seek medical advice.

Buffer TB contains chaotropes. It can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

PRODUCT DISCLAIMER

GeneAll® Exgene™ Soil DNA mini is for research use only, not for use in diagnostic procedure.

Product Specifications

| Specification | Exgene™ Soil DNA mini |
|---------------------------------------|-----------------------|
| Type | Spin |
| Maximum amount of starting samples | 500 mg soil sample |
| Maximum loading volume of spin column | 700 µl |
| Minimum elution volume | 30 µl |
| Maximum binding capacity | 100 µg |

Product Description

GeneAll® Exgene™ Soil DNA mini provides a convenient method for the isolation of total DNA from soil samples. This kit utilizes the powerful beads, the optimized buffer system and the advanced silica binding technology to purify nucleic acid suitable for many applications. These complex systems of this kit can deal with a number of different types of samples in the soil including plant tissues, bacteria, fungi spores and others. Also, it removes a humic acid and other PCR inhibitors from various soil samples efficiently. The humic acid, which is a sort of brownish colour, is a critical factor for soil treating experiments and if remained in eluate, this can have a negative effect on the DNA downstream applications.

GeneAll® Exgene™ Soil DNA mini provide a tube including powerful beads for strong pulverization. Soil samples are placed in this tube with lysis buffer, buffer SL, and crushed by bead-beater or vortex. After centrifugation, supernatant is mixed with precipitation buffer, buffer RH and buffer PD, to precipitate humic acid and protein. Then, the separated DNA part, supernatant, blend into the binding buffer, buffer TB, and DNA is bound on the silica membrane through centrifugation. Following washing step with buffer NW, the bound DNA is eluted by buffer EB. Purified DNA can be directly applicable in conventional PCR, restriction analysis, electrophoresis, and any other downstream applications.

PROTOCOL FOR

Exgene™ Soil DNA mini

- 1. Add up to 500 mg of soil sample to a Powerbead™ tube.**
- 2. Add 550 ul of buffer SL to the tube.**
- 3. Homogenize the sample in the Precellys® 24 (Bertin, France) equipment for twice of 23 seconds at 6500 rpm.**

Alternatively, secure tubes horizontally on a flat-bed vortex pad with tape and vortex at maximum speed for 10 minutes.
- 4. Centrifuge at $\geq 10,000 \times g$ for 10 minutes at room temperature and carefully transfer the supernatant to a 1.5 ml tube (provided).**
- 5. Add 50 ul of buffer RH.**
- 6. Add 300 ul of buffer PD and mix well by vortexing.**
- 7. Centrifuge at $\geq 10,000 \times g$ for 5 minutes at room temperature and carefully transfer the supernatant to a 2 ml tube (provided).**

Small pellet containing humic acid, cell debris, and protein can be formed in the collection tube after centrifugation. Be careful not to disturb this pellet.
- 8. Add 900 ul of buffer TB and mix well by vortexing.**

If buffer TB precipitation, pre-heat in a 56°C water bath to dissolve completely.
- 9. Transfer up to 700 ul of the mixture to a mini spin column.**
- 10. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.**

Discard the pass-through and reinsert the mini spin column back into the same tube.

11. Repeat two more times step 9 ~ 10 using the remainder of the sample.

12. Add 500 ul of buffer NW to the mini spin column.

13. Centrifuge at $\geq 10,000 \times g$ for 30 seconds at room temperature.

Discard the pass-through and reinsert the mini spin column back into the same tube.

14. Centrifuge at maximum speed for 1 minute at room temperature to remove residual wash buffer.

Transfer the mini spin column to a new 1.5 ml tube (provided).

Residual ethanol may interfere with downstream reactions. Care must be taken at this step for eliminating the carryover of buffer NW.

15. Add 50 ul of buffer EB to the center of the membrane in the mini spin column.

Incubate for 1 minute at room temperature. Centrifuge at $\geq 10,000 \times g$ for 1 minute at room temperature.

Elution volume can be decreased to 30 ul for high concentration of DNA, but this will slightly decrease in overall DNA yield. If maximum recovery of DNA is preferred or the starting materials contain large amount of DNA, elution can be done in 200 ul of buffer EB.

Troubleshooting for Exgene™ Soil DNA mini

| Facts | Possible Causes | Suggestions |
|--|---|---|
| <p>Low or no recovery</p> | <p>Too much starting material</p> | <p>Too much starting material lead to inefficient lysis, followed by poor DNA yields. Reduce the amount of starting material.</p> |
| | <p>Insufficient Homogenization</p> | <p>Check the step 3 of protocol. Insufficient homogenization time and condition is related to low recovery yield.</p> |
| <p>Low efficiency of DNA amplification</p> | <p>Excess amonut of template DNA</p> | <p>An excess amount of template DNA will inhibit a PCR reaction. The template DNA is needed to dilute.</p> |
| <p>Eluate does not preform well in the downstream application</p> | <p>Residual ethanol remains in eluate</p> | <p>To remove any residual ethanol included in buffer NW from mini spin column membrane, centrifuge again for complete removal of ethanol.</p> |
| <p>DNA eluate is brown</p> | <p>Humic acid is not be removed completely</p> | <p>With certain samples, a little humic acid can be remained in the eluate. In this case, we recommend using a GeneAll Expin Clean up SV kit to purify contaminated eluate.</p> |

Ordering Information

GeneAll® Hybrid-Q™ for rapid preparation of plasmid DNA

| Products | Type | Size | Cat. No. |
|-------------------|-------------|------|----------|
| Plasmid Rapidprep | mini / spin | 50 | 100-150 |
| | | 200 | 100-102 |

GeneAll® Exprep™ for preparation of plasmid DNA

| Products | Type | Size | Cat. No. |
|-------------------|---------------|-------|----------|
| Plasmid SV mini | spin / vacuum | 50 | 101-150 |
| | | 200 | 101-102 |
| | | 1,000 | 101-111 |
| Plasmid SV Midi** | spin / vacuum | 26 | 101-226 |
| | | 50 | 101-250 |
| | | 100 | 101-201 |

GeneAll® Exfection™ for preparation of highly pure plasmid DNA

| Products | Type | Size | Cat. No. |
|--------------------------------------|---------------|------|----------|
| Plasmid LE mini (Low Endotoxin) | spin / vacuum | 50 | 111-150 |
| | | 200 | 111-102 |
| Plasmid LE Midi* (Low Endotoxin) | spin / vacuum | 26 | 111-226 |
| | | 100 | 111-201 |
| Plasmid EF Midi* (Endotoxin Free) | spin | 20 | 121-220 |
| | | 100 | 121-201 |

GeneAll® Expin™ for purification of fragment DNA

| Products | Type | Size | Cat. No. |
|------------|----------------------|------|----------|
| Gel SV | mini / spin / vacuum | 50 | 102-150 |
| | | 200 | 102-102 |
| PCR SV | mini / spin / vacuum | 50 | 103-150 |
| | | 200 | 103-102 |
| CleanUp SV | mini / spin / vacuum | 50 | 113-150 |
| | | 200 | 113-102 |
| Combo GP | mini / spin / vacuum | 50 | 112-150 |
| | | 200 | 112-102 |

GeneAll® Exgene™ for isolation of total DNA

| Products | Type | Size | Cat. No. |
|------------------------|---------------|------|----------|
| Tissue SV mini* | spin / vacuum | 100 | 104-101 |
| | | 250 | 104-152 |
| Tissue SV Midi** | spin / vacuum | 26 | 104-226 |
| | | 100 | 104-201 |
| Tissue SV MAXI** | spin / vacuum | 10 | 104-310 |
| | | 26 | 104-326 |
| Tissue plus! SV mini* | spin / vacuum | 100 | 109-101 |
| | | 250 | 109-152 |
| Tissue plus! SV Midi** | spin / vacuum | 26 | 109-226 |
| | | 100 | 109-201 |
| Tissue plus! SV MAXI** | spin / vacuum | 10 | 109-310 |
| | | 26 | 109-326 |

GeneAll® Exgene™ for isolation of total DNA

| Products | Type | Size | Cat. No. |
|-------------------|---------------|------|----------|
| Blood SV mini | spin / vacuum | 100 | 105-101 |
| | | 250 | 105-152 |
| Blood SV Midi** | spin / vacuum | 26 | 105-226 |
| | | 100 | 105-201 |
| Blood SV MAXI** | spin / vacuum | 10 | 105-310 |
| | | 26 | 105-326 |
| Cell SV mini | spin / vacuum | 100 | 106-101 |
| | | 250 | 106-152 |
| Cell SV MAXI** | spin / vacuum | 10 | 106-310 |
| | | 26 | 106-326 |
| Clinic SV mini | spin / vacuum | 100 | 108-101 |
| | | 250 | 108-152 |
| Clinic SV Midi | spin / vacuum | 26 | 108-226 |
| | | 100 | 108-201 |
| Clinic SV MAXI** | spin / vacuum | 10 | 108-310 |
| | | 26 | 108-326 |
| Genomic DNA micro | spin | 50 | 118-050 |
| Plant SV mini | spin / vacuum | 100 | 117-101 |
| | | 250 | 117-152 |
| Plant SV Midi** | spin / vacuum | 26 | 117-226 |
| | | 100 | 117-201 |
| Plant SV MAXI** | spin / vacuum | 10 | 117-310 |
| | | 26 | 117-326 |
| GMO SV mini | spin / vacuum | 50 | 107-150 |
| | | 200 | 107-102 |
| Soil mini | spin | 50 | 114-150 |

GeneAll® GenEx™ for isolation of total DNA

| Products | Type | Size | Cat. No. |
|----------|----------------------------|------|----------|
| GenEx™ B | Sx [†] / solution | 100 | 220-101 |
| | | 500 | 220-105 |
| | | 100 | 220-301 |
| GenEx™ C | Sx [†] / solution | 100 | 221-101 |
| | | 500 | 221-105 |
| | | 100 | 221-301 |
| GenEx™ T | Sx [†] / solution | 100 | 222-101 |
| | | 500 | 222-105 |
| | | 100 | 222-301 |

GeneAll® DirEx™ Single tube DNA extraction buffer for PCR

| Products | Type | Size | Cat. No. |
|----------|----------|------|----------|
| DirEx™ | solution | 50 | 250-050 |

| Products | Type | Size | Cat. No. |
|---|-------------------|---------|----------|
| GeneAll® RNA Series for preparation of RNA | | | |
| RiboEx™ | solution | 100 | 301-001 |
| | | 200 | 301-002 |
| Hybrid-R™ | spin | 100 | 305-101 |
| Hybrid-R™ Blood RNA | spin | 50 | 315-150 |
| Hybrid-R™ miRNA | spin | 50 | 325-150 |
| RiboEx™ LS | solution | 100 | 302-001 |
| | | 200 | 302-002 |
| Riboclear™ | spin | 50 | 303-150 |
| Ribospin™ | spin | 50 | 304-150 |
| Ribospin™ vRD | spin | 50 | 302-150 |
| Ribospin™ Plant | spin | 50 | 307-150 |
| Allspin™ | spin | 50 | 306-150 |
| GeneAll® AmpONE™ for PCR amplification | | | |
| Taq DNA polymerase | (2.5 U/μℓ) | 250 U | 501-025 |
| | | 500 U | 501-050 |
| | | 1,000 U | 501-100 |
| α-Taq DNA polymerase | (2.5 U/μℓ) | 250 U | 502-025 |
| | | 500 U | 502-050 |
| | | 1,000 U | 502-100 |
| Pfu DNA polymerase | (2.5 U/μℓ) | 250 U | 503-025 |
| | | 500 U | 503-050 |
| | | 1,000 U | 503-100 |
| Hotstart Taq DNA polymerase | (2.5 U/μℓ) | 250 U | 531-025 |
| | | 500 U | 531-050 |
| | | 1,000 U | 531-100 |
| Clean Taq DNA polymerase | (2.5 U/μℓ) | 250 U | 551-025 |
| | | 500 U | 551-050 |
| | | 1,000 U | 551-100 |
| Clean α-Taq DNA polymerase | (2.5 U/μℓ) | 250 U | 552-025 |
| | | 500 U | 552-050 |
| | | 1,000 U | 552-100 |
| Taq Master mix | 0.5 ml x 2 tubes | 2x | 511-010 |
| | 0.5 ml x 10 tubes | 2x | 511-050 |
| α-Taq Master mix | 0.5 ml x 2 tubes | 2x | 512-010 |
| | 0.5 ml x 10 tubes | 2x | 512-050 |

* GeneAll® Tissue SV mini, Midi, and MAXI plus! kit provide the additional methods for the purification from animal whole blood.

** GeneAll® SV Midi / MAXI kits require the centrifuge which has a swinging-bucket rotor and ability of 4,000 ~ 5,000 xg.

† On the basis of DNA purification from 300 ul whole blood, 2 x 10⁶ cells or 10 mg animal tissue.

†† On the basis of DNA purification from 10 ml whole blood. 1 x 10⁶ cells or 100 mg animal tissue.

| Products | Type | Size | Cat. No. |
|--|-------------|----------------|----------|
| GeneAll® AmpONE™ for PCR amplification | | | |
| Taq Premix | 96 tubes | 20 μℓ | 521-200 |
| | | 50 μℓ | 521-500 |
| α-Taq Premix | 96 tubes | 20 μℓ | 522-200 |
| | | 50 μℓ | 522-500 |
| Taq Premix (w/o dye) | 96 tubes | 20 μℓ | 524-200 |
| α-Taq Premix (w/o dye) | 96 tubes | 20 μℓ | 525-200 |
| dNTP mix | 2.5 mM each | 500 μℓ | 509-020 |
| dNTP set (set of dATP, dCTP, dGTP and dTTP) | 100 mM | 1 ml x 4 tubes | 509-040 |

* Each dNTP is available



GeneAll

GENEALL BIOTECHNOLOGY CO., LTD

www.geneall.com

GeneAll Bldg., 128 Oguem-dong,
Songpa-gu, Seoul, KOREA 138-859

E-MAIL sales@geneall.com

T E L 82-2-407-0096

F A X 82-2-407-0779

©2012 GeneAll Biotechnology, All right reserved

Edited by SR
Designed by Park Eun Ah